



PAT LAL - True Low Air Loss Support Surface

Therapeutic Support Surface prevents and treats pressure ulcers by molding to the body to maximize contact, redistributing weight as uniformly as possible, and reducing pressure to below capillary closure (32mm hg)



"Up to 1000 liters of air flow per minute delivered to the surface provides an ideal microclimate for wound healing as the surface gently pulsates to aid in stimulating capillary blood flow."

Features

- ✓ True Low Air loss therapy provides the ideal microclimate for wound healing
- ✓ Three pressure zones provide advanced pressure relief
- ✓ Gentle pulsation to assist in stimulating blood flow
- ✓ Comfortable low shear, low friction, vapor permeable cover
- ✓ One touch AUTOFIRM to facilitate transfers and position changes
- ✓ Keypad lockout eliminates unauthorized changes to therapeutic settings
- ✓ All Telemade air systems are available in bariatric widths and weight capacities.



Nu Dimension Tech Ventures LLC 3215 Huffman Eastgate Rd Huffman TX 77336

Website: www.telemade.com Phone: 281-324-2900 Email: mayre@tele-made.com