




PAT-LAL

True Low Air Loss Support Surface



- True low air loss therapy provides the ideal microclimate for wound healing
- Three pressure zones provide advanced pressure relief
- Gentle pulsation to assist in stimulating blood flow
- Comfortable low shear, low friction, vapor permeable cover
- One touch AUTOFIRM to facilitate transfers and position changes
- Key pad lockout eliminates unauthorized changes to therapeutic settings



“Up to 1000 liters of air flow per minute are delivered to the patient surface to provide an ideal microclimate for wound healing as the surface gently pulsates to aid in stimulating capillary blood flow.”

The Telemade, Patriot LAL provides advanced pressure relief with true low air loss therapy to provide the ideal microclimate for wound healing. Delivering up to 1000 liters of air per minute to the true low air loss surface reduces heat and moisture build-up at the patient, surface interface. Reducing heat and moisture reduces the detrimental effects they can have on pressure, shear and friction. Three modes of operation and simple comfort adjustments allow the user to customize therapy. Low Air Loss mode provides continuous advanced pressure relief with true low air loss. Pulsation mode provides all the benefits of the Low Air Loss mode with the addition of gentle pulsation to aid in stimulating capillary blood flow. By engaging the Autofirm mode the Patriot LAL quickly provides a uniform firm surface to facilitate care giving from position changes to transfers.



www.tele-made.com 281-324-2900