




PAT-ALAL

True Low Air Loss Therapy with Alternating Pressure



- True low air loss therapy provides the ideal microclimate for wound healing
- Three pressure zones provide advanced pressure relief
- Gentle alternation to assist in stimulating capillary circulation
- Comfortable low shear, low friction, vapor permeable cover
- One touch AUTOFIRM to facilitate transfers and position changes
- Key pad lockout eliminates unauthorized changes to therapeutic settings



"Up to 1000 liters of air flow per minute are delivered to the patient surface to provide an ideal microclimate for wound healing."



Gentle alternation enhances capillary circulation while maintaining constant advanced pressure relief

The Telemade, Patriot ALAL provides advanced pressure relief with true low air loss therapy to provide the ideal microclimate for wound healing. Delivering up to 1000 liters of air per minute to the true low air loss surface reduces heat and moisture build-up at the patient, surface interface. Reducing heat and moisture reduces the detrimental effects they can have on pressure, shear and friction. Three modes of operation and simple comfort adjustments allow the user to customize therapy. Static mode provides continuous advanced pressure relief with true low air loss. Alternation mode provides all the benefits of the Static mode with the addition of gentle alternating pressure to aid in stimulating blood flow. By engaging the Autofirm mode the Patriot ALAL quickly provides a uniform firm surface to facilitate care giving from position changes to transfers.

PAT-ALAL HCPCS: E0277

Manufactured by Telemade
3215 Huffman Eastgate, Huffman, TX 77336

Telemade[®]
Specialty Medical Products
Made in the U.S.A.
(281) 324-2900