




PAT-LAL

Cost-effective, True Low Air Loss Therapy with Pulsation



- True low air loss therapy provides the ideal microclimate for wound healing
- Three pressure zones provide advanced pressure relief
- Gentle pulsation to assist in stimulating blood flow
- Comfortable low shear, low friction, vapor permeable cover
- One touch AUTOFIRM to facilitate transfers and position changes
- Key pad lockout eliminates unauthorized changes to therapeutic settings



“Up to 1000 liters of air flow per minute are delivered to the patient surface to provide an ideal microclimate for wound healing as the surface gently pulsates to aid in stimulating capillary blood flow.”

The Telemade, Patriot LAL provides advanced pressure relief with true low air loss therapy to provide the ideal microclimate for wound healing. Delivering up to 1000 liters of air per minute to the true low air loss surface reduces heat and moisture build-up at the patient, surface interface. Reducing heat and moisture reduces the detrimental effects they can have on pressure, shear and friction. Three modes of operation and simple comfort adjustments allow the user to customize therapy. Static mode provides continuous advanced pressure relief with true low air loss. Pulsation mode provides all the benefits of the Static mode with the addition of gentle pulsation to aid in stimulating blood flow. By engaging the Autofirm mode the Patriot LAL quickly provides a uniform firm surface to facilitate care giving from position changes to transfers.

Specifications:	Mattress	Control Unit
	Dimensions: 9" x 36" x 80"	Static, Pulse and AutoFirm Operating Modes
	8" cushions over 1" foam base	Up to 1200 liters per minute air flow
	Comforter: Air diffusing microfiber urethane coated nylon taffeta	120 V, 60Hz
	Base: Heavy duty vinyl, 10 security straps	
	Patient Weight Limit 350 lb.	

Manufactured by Telemade, Huffman, Texas 77336 (281) 324-2900



Specialty Medical Products
Made in the U.S.A.
(281) 324-2900