



## Low Air Loss - Alternating Pressure Relief Therapy


### LAL-AP

The Patriot LAL-AP delivers proven alternating pressure relief therapy with low air loss. Two internal pumps work in unison to provide continuous smooth alternation to the patient surface as excess air is redirected to the vented cover sheet to aid in reducing heat and moisture build-up. Internal cushion pressures are continually monitored by the microprocessor to ensure proper inflation and therapeutic alternating pressure relief.




## Therapy Support Surface

- 8" Alternating pressure therapy surface
- 20 Individual nylon, urethane reinforced air cells
- Advanced Velcro and cushion sleeve securely keeps cushions in place.
- Dedicated hose directs air loss to the vented comforter
- Low shear, low friction nylon taffeta patient surface
- Comforter is water-proof yet vapor permeable
- Stays inflated in temporary power outages and has internal foam base for added power loss protection



## Control Unit



- Simple user friendly controls for ease of use
- Soft touch control panel
- Internal monitor continually monitors and maintains preset pressures
- Auto Firm, Low Air Loss and Alternating Modes
- Key pad lock-out
- Will stay inflated during temporary power outages

Mattress: 80" x 36" x 8"

Max Weight Limit: 350 lb.

  
**281-324-2900**