




## Low Air Loss/Alternating Pressure Mattress System

The Patriot LAL-AP combines dual therapies of low air loss and alternating pressure in one unit. Ideal for patients who require treatment for pressure ulcers stages I-IV, the Patriot Digital pump alternates in 10 minute intervals and the adjustable LED digital air flow touch panel provides custom pressure settings for each patients individual requirements. (Medicare code E0277)



### Patriot PAT-LAL-AP Control Unit

- The Patriot PAT-LAL-AP digital controller provides a simple user friendly interface for full mattress operation and proven alternating low air loss therapy.
- Alternating Pressure-Low Air Loss.
- Digital adjustable LED weight setting allows pressure to be customized according to each patient's therapeutic requirements.
- System continually monitors clinicians setting and automatically adjusts to maintain proper pressure.
- Static function for safe and easy patient transfers and position changes.
- Static Mode stops Alternation and provides low air loss pressure redistribution.
- Audio/visual low pressure alarm alerts clinicians of potential compromised airflow.
- Model: PAT-LAL-AP
- 20, alternating air cells provide maximum pressure redistribution.
- Air cells are securely fastened to the base via unique hook and loop system as well as a center cushion coupler to eliminate cushion migration.
- The Patriot LAL-AP low shear/low friction nylon taffeta comforter cover is fluid resistant, anti-microbial/anti-bacterial and vapor permeable.
- Integrated high density foam base provides protection in case of unexpected power outage.
- Weight capacity 400 lbs.
- 2 year warranty.
- Made in the USA.
- 36" x 80" x 8". Custom sizes available.

